

Medical Trolley

JMT500

Specifications



Power System	220V AC input, 12V DC output
Power status indicator	LED green-orange-red
Monitor arm	Gas spring
Monitor weight	5-10kg
VESA mount	75x75mm & 100x100mm
Lift mechanism	Motorized (Push 500N)
Height adjustable	Desktop 80cm to 110cm Monitor 112cm to 153cm
Main Material	Aluminium alloy/Stainless steel/MDF with medical grade coating
Mounting option	DIN rails left and right
Work surface	510mm (L) x 580mm(W)
Keyboard tray	260mm (D) x 536mm(W)
Transport handle	Front
Wheels	4 castors, 2 lockable castors in front castor dimension: 125mm
Battery	LiFePO 40Ah
Battery compartment	180mm(L) x 290mm(W) x 159mm(H)
Pre-assembled cables	3 x 2.0 USB cables 1 x LAN cable (RJ45 8P8C M to F)
Power cable	3 meter retractable power cable (220-240V/50-60Hz/6-8A)
Compliance	CE, RoHS, NEN/EN/IEC 3140, NEN/EN/IEC 62353, Class 2 double insulation
Warranty	Mechanical parts: 2 years Power system/electrical parts: 2 year
<div><div>MD</div><div>This is a medical device as defined in MDR 2017/745</div></div>	

Fusion Medical

Web: www.fusionmedical.co.za

Email: info@fusionmedical.co.za

Tel: +27 82 469 8816

Address: 20 Felicia Street, Brackenhurst,
Alberton, 1448, Gauteng, South Africa



FUSION MEDICAL
TECHNOLOGIES



娅诺
Jannu